FR

COMPOSITION

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INDICATION

Permanent root canal filling and sealing in dental procedures, with systematic insertion of points.

PROPERTIES

To act as perfect root canal filling materials, pastes should be neither resorbable once hardened nor retractable over the years.

Their antiseptic and anti-inflammatory activity should last for several hours after being placed in the canal.

These pastes should be radiopaque and easy to insert in the canal. ENDOMETHASONE N Powder has all the various above-described qualities, owing to its excipients and its antiseptic active substance (thymol iodide).

DIRECTIONS FOR USE

The treated tooth should be perfectly isolated (using a cotton roll or a rubber dam).

After carefully shaping the canal, cleaning it with an antiseptic solution, rinsing it and drying it with paper points, prepare the canal paste extemporaneously, progressively pouring the powder into the liquid (eugenol or ENDOMETHASONE Liquid) at a ratio of 2 spoonfuls of powder to 3-6 drops of liquid, until the paste has the consistency required for the selected filling technique:

- using a paste carrier

- and/or by coating cones with a cement sealer for condensation techniques.

Then fill the root canal using a paste-carrier and position into the paste one or several points, which have previously been adapted to the canal shape.

Take X-rays to ensure that root canals have been thoroughly filled.

INTERACTIONS

Do not use in combination with composite materials, since eugenol-based preparations inhibit their polymerisation.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

As with any other root canal sealer, a canal can only be sealed if it has previously been reamed, cleaned and disinfected thoroughly.

Determine the working length of the root canal before filling so as to prevent any apical extrusion, since the paste is not resorbable.

Risk of allergy to any of the constituents, particularly in case of contact with the gums. In case of accidental contact, rinse with plenty of water.

Due to electrostatics, ENDOMETHASONE, as all powders, tends to agglomerate in the bottle. However, the capacity mentioned on the label corresponds to the weighed contents. It is therefore important to shake the bottle prior to using the measuring spoon, in order to loosen the powder.

C - Corrosive (contains hydrocortisone acetate and thymol diiodide).
N - Dangerous for the environment.

Danger of cumulative effects.

Causes burns.

Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable protective clothing, gloves and eye/face protection.

In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Avoid release to the environment. Refer to special instructions/Safety data sheets.

This material and its container must be disposed of as hazardous waste.

STORAGE

Do not store above 25°C and protect from moisture.

PRESENTATIONS

14 g or 42 g bottle of powder + 1 measuring spoon.

3 x 14 g bottles of powder + 1 measuring spoon.

Kit : 1 x 14 g bottle of powder + 1 x 10 ml bottle of liquid + 1 measuring spoon.

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Eugenol 91.0 ml
Excipient q.s.ad. 100.0 ml

INDICATION

Liquid to be mixed with Endomethasone Powder to form a root-canal filling paste.

INSTRUCTIONS FOR USE

Prepare the root canal filling paste extemporaneously, progressively mixing 2 spoonfuls of ENDOMETHASONE Powder with 3 to 6 drops of ENDOMETHASONE Liquid,

until the paste reaches the consistency required for the chosen root-canal filling technique.

Fill the root canals using a paste-carrier and insert into the paste one or several points after adapting it/them to the canal shape.
Take X-rays to ensure that the canals have been thoroughly filled.

WARNINGS AND PRECAUTIONS FOR USE

Eugenol inhibits the polymerization of resin materials.

Xi – Irritating (contains eugenol).
Irritating to skin. May cause sensitization by skin contact.
Avoid contact with skin. Wear suitable gloves.

PRESENTATION

10 ml bottle.

For professional use only in the practice of dentistry.

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