F

COMPOSITION

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UK

COMPOSITION

Acetone	37.60 g
Ethyl acetate	42.80 g
Excipients q.s.	100.00 ml

INDICATIONS

Alcohol or chloroform which have been used so far to dry dental cavities, have two main disadvantages:

1. They do not always remove the grease present in cavities.

2. They are contraindicated when using some filling or sealing materials (porcelains, resins, etc.).

HYDROL simultaneously removes moisture and all traces of grease.

HYDROL can be used with any type of filling or sealing material.

HYDROL is economical in use and permanently solves one of the practitioner's constant problems.

DIRECTIONS FOR USE

Dab the area to be dried with a cotton ball soaked in HYDROL. Wait a few seconds for complete evaporation of the liquid.

SPECIAL WARNINGS

F - Highly flammable

IMPORTANT REMARKS

Before filling the tooth with a temporary self-curing paste or a permanent acrylic resin filling material, wait until the cavity is completely dry.

HYDROL is not designed for degreasing the gingival mucosa, as prolonged application may cause a burning sensation.

Carefully close the bottle after use.

PRESENTATIONS

45 ml bottle.
250 ml bottle.

For professional dental use only.

Distributor: SEPTODONT Ltd.

Unit R & S - Orchard Business Centre
St Barnabas Close
ALLINGTON, Maidstone, KENT ME16 OJZ - UK

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