F

## PROPRIÉTÉS

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UK

## PROPERTIES

PROVISCCELL, an eugenol-free temporary cement is simple to apply and can be easily removed from plastic and metal crowns. PROVISCCELL has no influence on the hardening of resin-based sealing or filling materials.

## INDICATION

Temporary sealing.

## DIRECTIONS FOR USE

### • Tubes

PROVISCCELL should be used at room temperature (18 - 23 °C) and at a normal humidity level (45-55 %). In case a cool glass slab is used, the working time will be longer.

- Squeeze equal amounts of base and catalyst into a mixing pad or glass plate.
- Mix the two pastes for approximately 20 seconds to a homogeneous consistency. Use a plastic or rustproof metal spatula made for mixing.
- Apply a thin layer of PROVISCCELL inside the temporary crown or bridge and cement it to the prepared tooth. PROVISCCELL flows freely while the crown or bridge is being put in place.
- PROVISCCELL remains workable for about 1 1/2-2 minutes and sets in the mouth within 4-6 minutes.
- 3-5 minutes after PROVISCCELL has set, remove excess material.

### Remarks

To obtain a softer material and a longer setting time add some more base paste or a small amount of petrodatum.

A harder and faster setting material is obtained by adding more catalyst paste. A moist environment (humid stumps) facilitates the removal of temporary cement.

### • Syringes

- Remove the cap by turning 90° anticlockwise and discard it.
- Attach a mixing tip to the syringe. Line up the tip notch and the syringe notch.
- Turn the mixing tip 90° clockwise until it locks.
- The cement is automatically mixed in the tip and can be applied directly.
- At first use, extrude pea-sized amount of material and discard it.
- Apply the mixed cement into the dried temporary restoration.
- Fix the temporary restoration onto the dried prepared core, exerting gentle pressure.
- After setting (approximately 4 minutes), remove excess cement with a probe.

**Leave the used mixing tip on the syringe to serve as a cap until next use.**

## SPECIAL WARNINGS

In case of exposed pulp or very low dentin thickness, it is recommended to use a pulp cap before. Do not use PROVISCCELL if the patient is known to be allergic to fatty acids and colophony.

## STORAGE

Keep the tubes at a temperature below 23 °C and close them well after use.

Keep the syringes at a temperature below 25°C.

## PRESENTATIONS

### • Tubes

Set containing :

One 25 g tube of base,

One 25 g tube of catalyst,

One mixing pad.

### • Syringes

2 x 11 g syringes + 20 mixing tips

50 mixing tips.

*For professional use only in the practice of dentistry.*

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