## *Resorcell Gaze*

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Knitted oxidized cellulose gauze (pure cotton).

## PROPERTIES

It is a white, oxidized cellulose gauze, which can be cut without fraying.

After the gauze has been saturated with blood, it swells into a red, gelatinous mass, which aids in the formation of a clot, thereby serving as a haemostatic adjunct in the control of haemorrhage.

When used properly, in minimal amounts, the gauze is absorbed in the implantation site.

Its bacteriostatic/bactericidal properties inhibit aerobic and anaerobic bacteria cell growth and multiplication. However, this product is not intended as a substitute for systemically administered antibiotics, to fight or prevent postoperative infections.

## INDICATIONS

Haemostatic adjunct in the control of capillary and small vascular bleeding, after tooth extraction.

## CONTRAINDICATIONS

Do not use for implantation in bone tissue when irrigation with blood is not sufficient, as the acidity of the product causes a decrease in callus formation.

## INSTRUCTIONS FOR USE

Sterile technique should be observed when removing the product from its packaging.

Use only as much gauze as is necessary for haemostasis, holding it firmly in place until bleeding stops. Absorption depends upon several factors including the amount use, degree of saturation with blood and nature of the operative site. Furthermore, the gauze should preferably be packed loosely in the socket, to prevent the risk of tissue necrosis due to extra pressure.

## SPECIAL WARNINGS AND PRECAUTIONS FOR USE

Remove any excess gauze after haemostasis and prior to surgical closure to facilitate absorption and minimize the risk of encapsulated cyst formation.

The gauze should be applied dry and should not be moistened prior to application. It should not be impregnated with anti-infective agents or with any other materials such as buffering or haemostatic substances. Its haemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.

Check package integrity before use.

 **Do not use if packaging has been damaged or opened before first use.**

Do not use gauze which may have become contaminated, since the product cannot be re-sterilized.

## ADVERSE REACTIONS

May cause burning or stinging sensations due to the low pH of the product.

Do not use as a dressing for superficial wounds, as this product inhibits epithelization.

## STORAGE

Protect from direct sunlight and keep in the original packaging.

Store dry between 0°C and 10°C.

Do not use beyond the expiry date stated on the label. Browning of the gauze and tendency to crumble indicate improper storage.

## PRESENTATION

Carton containing 10 gauzes (5 x 1.25 cm).

Gauze in a double peel-pack, sterilized by gamma irradiation.

**For professional use only in the practice of dentistry.**

Distributor : SEPTODONT Ltd.

Units R & S – Orchard Business Center  
St Barnabas Close Allington, Maidstone  
Kent ME16 0JZ - UK

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