Biodentine™ XP 500

Biodentine™ XP 200

Active Biosilicate Technology™

en	Bioactive dentine substitute	1
fr	Substitut dentinaire bioactif	4
es	Sustituto bioactivo de la dentina	7

Biodentine™ XP is a bioactive dentine substitute from the “Active Biosilicate Technology™” innovation.

COMPOSITION

Powder part is composed of tricalcium silicate, zirconium oxide, calcium oxide, calcium carbonate and iron oxides.

Liquid part is composed of calcium chloride, polycarboxylate and purified water

INTENDED USE

Intended purpose

Biodentine™ XP is a bioactive dentine substitute.

Clinical benefits

Preservation of the tooth vitality (except for apexification, perforations and furcations);

Preservation of the tooth functionality;

Intended population

Biodentine™ XP is intended to be applied on temporary teeth (from 2 years old), and on permanent teeth (mature or immature) depending on indication. Please refer to the dedicated section (§Indications for use).

Intended users

Dental healthcare professional only.

Indications for use

In the crown:

For temporary teeth (children from 2 years of age) and permanent teeth (immature or mature):

- Permanent dentine restoration under composites or Inlay/Onlay;
- Temporary dentine-enamel restoration;
- Restoration of deep and/or large coronal carious lesions (sandwich technique);
- Restoration of cervical or radicular lesions.

On the pulp

For temporary teeth (children from 2 years of age) and permanent teeth (immature or mature):

- Pulp capping (direct and indirect);
- Pulpotomy for diagnosed symptoms of reversible pulpitis and irreversible pulpitis where bleeding is controlled within 5 minutes.

In the root

For permanent teeth (immature or mature):

- Repair of root perforations;
- Repair of furcation perforations;
- Repair of perforating internal resorptions;
- Repair of external resorption;
- Root-end filling in endodontic surgery (retrograde filling).

For immature permanent teeth

- Apexification;

CONTRA INDICATION

- Hypersensitivity to one of the ingredients.
- Restoration of large loss of tooth substance subjected to high stresses.

PERFORMANCE CHARACTERISTICS

Technical specifications

- Compressive strength: $\geq 150\text{MPa}$
- Working time: $\geq 1\text{ min}$
- Setting time: the average setting time of the product is 12 minutes, starting from the end of mix. Depending on some conditions, that include storage time and conditions, the setting time may vary (from 9 to 25 minutes).

Technical Performances

Biodentine™ XP is characterized by:

- its sealing ability by creating a tight seal within the dentin tubules;
- its bioactivity with hydroxyapatite formation;
- its dentine substitution role.

Clinical performance

- Preservation of the pulp;
- Preservation of the amelo-dentinal structure;
- Preservation of the root canal structure;
- Root edification;

EXPECTED LIFE-TIME

Biodentine™ XP is in contact with pulpal, dentinal, or enamel tissues during 10 years.

INSTRUCTION FOR USE

CAUTION: All the following instructions must be carried out within the working area.

1. Take out the cartridge from the protection foil.
2. Insert the cartridge into the Biodentine™ Mixer activation socket and lock the cartridge making a 90° clockwise turn.
3. Activate the cartridge using Biodentine™ Mixer. Low down the lever completely until you hear the clicking sound.
4. Unlock the cartridge from the activation socket making an anti-clockwise 90° turn.
5. Immediately after the activation step, move the cartridge up and lock it on the Biodentine™ Mixer cartridge holder making a clockwise 90° turn.
6. Close the lid and push on the “Start” button.
7. At the end of mix, open the lid, unlock the cartridge from its holder making an anti-clockwise 90° turn.
8. Immediately after mixing hold Biodentine™ Gun, raise the lever until it stops in horizontal position, align the cartridge with the Biodentine™ Gun notches. Insert the cartridge into the notches and lock the cartridge making a clockwise 90° turn.
9. Once the cartridge locked in the gun, the cartridge nozzle can be bent manually in order to make the placement of the material easier.

CAUTION: Do not bend the nozzle more than 4 times at 90° to avoid any breakage of the nozzle.

10. Hold Biodentine™ Gun with both hands, push down the lever until the material appears in the nozzle.
11. Hold Biodentine™ Gun with one hand, pushing down on the lever to allow the cement to flow out of the cartridge. Apply the cement directly into to the affected tooth cavity with the Biodentine™ Gun or with an adequate dental instrument.
12. After use, raise the lever until it stops in horizontal position, remove the cartridge from Biodentine™ Gun making an anti-clockwise 90° turn and throw it in the regular bin.

WARNINGS

- If the cartridge falls on the floor, clean it with ethanol to avoid contamination.
- The device may cause eye irritation and skin burn Do not manipulate the device without appropriate personal protective equipment.
- Use a rubber dam to ensure that the cement is not in contact with the mucosa or swallowed.
- Before performing the permanent enamel restoration or removing the matrix (when used); check that the product is hard enough.
- Single-use product: do not use the mixed paste on several patients to avoid cross contamination.
- Biodentine™ XP is only intended to be mixed with Biodentine™ Mixer and applied with Biodentine™ Gun. Otherwise, technical specifications are not guaranteed.
- Flush with sterile saline solution to remove any residue of irrigants and dry without totally dehydrating the root canal for endodontic procedures.

SIDE-EFFECTS

There are no known side effects.

DISPOSAL

No specific recommendations to dispose.

STORAGE

Store between 5°C - 25°C (41°F to 77°F) and protect from moisture.

CONTENT

Box containing 10 cartridges of 700 mg of powder and 190 µL of liquid.

Box containing 10 cartridges of 468 mg of powder and 133 µL of liquid.

	Storage temperature range
	Consult instructions for use
	Caution
	Keep dry, protect from excessive humidity, rain or damp conditions
	Do not reuse
Rx only	By prescription only
	Catalogue or model number
	Manufactured by
	Distributed by

CAUTION: Federal law restricts this device to sale by or on order of a dentist

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SEPTODONT Inc