Pediatrie

**ACTIVE
BIOSILICATE
TECHNOLOGY**

BiodentineTM

Dentinesubstituut

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Authors	Title	Journal	Year	Ref.
Kurun Aksoy M, Tulga Oz F, Orhan K.	Evaluation of calcium (Ca ²⁺) and hydroxide (OH ⁻) ion diffusion rates of indirect pulp capping materials	International Journal of Artificial Organs	2017	1
Gong V, França R.	Nanoscale chemical surface characterization of four different types of dental pulp-capping materials	Journal of Dentistry	2017	3
Costa F, Sousa Gomes P, Fernandes MH.	Osteogenic and Angiogenic Response to Calcium Silicate-based Endodontic Sealers	Journal of Endodontics	2016	4
El Meligy OA, Allazzam S, Alamoudi NM.	Comparison between Biodentine™ and formocresol for pulpotomy of primary teeth: A randomized clinical trial	Quintessence	2016	
Grewal N, Salhan R, Kaur N, Patel HB.	Comparative evaluation of calcium silicate- based dentin substitute (Biodentine™) and calcium hydroxide (pulpdent) in the formation of reactive dentin bridge in regenerative pulpotomy of vital primary teeth: Triple blind, randomized clinical trial	Contemporary Clinical Dentistry	2016	
Özyürek T, Demiryürek EÖ.	Comparison of the antimicrobial activity of direct pulp-capping materials: Mineral trioxide aggregate-Angelus and Biodentine™	Journal of Conservative Dentistry	2016	8
Martens L, Rajasekharan S, Cauwels R.	Endodontic treatment of trauma-induced necrotic immature teeth using a tricalcium silicate-based bioactive cement. A report of 3 cases with 24-month follow-up.	European Journal of Pediatric Dentistry	2016	
Evren OK, Altunsoy M, Tanriver M, Capar ID, Kalkan A, Gok T.	Fracture resistance of simulated immature teeth after apexification with calcium silicate-based materials	European Journal of Pediatric Dentistry	2016	
Nowicka A, Wilk G, Lipski M, Kolecki J, Buczkowska-Radlińska J.	Tomographic Evaluation of Reparative Dentin Formation after Direct Pulp Capping with Ca(OH) ₂ , MTA, Biodentine™, and Dentin Bonding System in Human Teeth	Journal of Endodontics	2015	2
Camilleri J.	Staining Potential of Neo MTA Plus, MTA Plus, and Biodentine™ Used for Pulpotomy Procedures	Journal of Endodontics	2015	7
Koubi G, Colon P, Franquin JC, Hartmann A, Richard G, Faure MO, Lambert G.	Clinical evaluation of the performance and safety of a new dentin substitute, Biodentine™, in the restoration of posterior teeth - a prospective study	Clinical Oral Investigation	2012	6
Peng L1, Ye L, Tan H, Zhou X.	Evaluation of the formocresol versus mineral trioxide aggregate primary molar pulpotomy: a meta-analysis	Oral Surgery, Oral Medicine, Oral Pathology, Oral Radiology	2006	5

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