
GenENDO RevO-S+

Niti Root Canal Instruments

- fr - Évasement coronaire et mise en forme canalaire
- en - Coronal flaring and root canal shaping
- es - Ensanchamiento coronal y conformación de los conductos radiculares
- de - Koronale Erweiterung und Wurzelkanalaufbereitung
- it - Svasatura coronale e sagomatura di canali radicolari
- nl - Verwijding van het coronale gedeelte en vormgeving van het wortelkanaal
- pl - Poszerzenie części koronowej i kształtowanie kanałów koronowych
- cs - Koronální flaring a preparace kořenového kanálku
- tr - Koronal genişletme (flaring) ve kök kanal şekillendirme
- ar - توسيع القنوات التاجية وتثبيتها

CE 0459

	O	%	L	Stop	Ring	RPM	Recommended Torque (N.cm)
RevO-S+ SC1	25	6	21	BLACK	WHITE	250 - 400	1.8 - 2.5
RevO-S+ SC2	25	4	21 - 25	GREY	YELLOW	250 - 400	1.8 - 2.5
RevO-S+ SU	25	8	21 - 25	3 BLACK	RED	250 - 400	1.8 - 2.5

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
- Caractéristiques**
Type 2 rotating instruments under EN ISO 3630-1:2019 standard:
Maximum recommended torque (not required): see table "Recommended Torque".
Material of the operative part: Nitinol.
Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (Type 1).
Number of uses: maximum of 5 canals recommended if the instrument is not visibly damaged.
Medical device class according to Directive 93/42/EEC and MDR 2017/745: IIa.
- Protocol**
When increased penetration of a manual file is possible up to apical constriction upon initial penetration, using RevO-S+ SC1 is unnecessary and RevO-S+ SC2 should be used directly.
RevO-S+ SC1 is used with a "wave" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters or until a constraint is felt. At each end of progression, a series of alternating partial support movements can be performed. The initial canal portions - cervical and middle - are thus successively shaped using a "crown-down" approach.
RevO-S+ SC2 and RevO-S+ SU are used with a so-called "wavy" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters. When the instrument reaches the apical limit of preparation, a series of selective accent brushing movements should be performed to finalize cleaning and shaping of the instrument.
- Operating procedure**
Where increased penetration of a manual file is possible up to apical constriction upon initial penetration, using RevO-S+ SC1 is unnecessary and RevO-S+ SC2 should be used directly.
RevO-S+ SC1 is used with a "wave" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters or until a constraint is felt. At each end of progression, a series of alternating partial support movements can be performed. The initial canal portions - cervical and middle - are thus successively shaped using a "crown-down" approach.
RevO-S+ SC2 and RevO-S+ SU are used with a so-called "wavy" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters. When the instrument reaches the apical limit of preparation, a series of selective accent brushing movements should be performed to finalize cleaning and shaping of the instrument.
- Protocol**
When increased penetration of a manual file is possible up to apical constriction upon initial penetration, using RevO-S+ SC1 is unnecessary and RevO-S+ SC2 should be used directly.
RevO-S+ SC1 is used with a "wave" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters or until a constraint is felt. At each end of progression, a series of alternating partial support movements can be performed. The initial canal portions - cervical and middle - are thus successively shaped using a "crown-down" approach.
RevO-S+ SC2 and RevO-S+ SU are used with a so-called "wavy" motion using three successive movements, from coronal to apical, centered in the canal, allowing progression of a few millimeters. When the instrument reaches the apical limit of preparation, a series of selective accent brushing movements should be performed to finalize cleaning and shaping of the instrument.

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- Inspection de reprocessing**
For all metal devices, the use of antiseptic disinfecting and cleaning agents is recommended.
For your own safety, please wear personal protective equipment (gloves, glasses and mask).
Do not use cleaning or disinfection agents containing phenol, aldehyde and alkaline composition.
Always follow the instructions for use provided by the manufacturer of the products.
- Storage**
Store the products in a dry, clean and dust-free environment at the temperature specified by the paper-plastic pouch by the steam sterilizer manufacturer.
Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete re-use should be performed.
- Disposal**
After use, instruments must be placed in a secure container, used for collecting cutting instruments (like needles or disposable bistouris) as per good dentistry practice.

- 11. Symbols**
NITI raw material
Quantity
CLASSICS' handle
Medical device

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