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Referenties

- 1 Bastos, M.C., Araujo Albuquerque, F. G., Cetira Felho, E. L. et al. (2024) 'Clinical and radiographic success of pulpotomy and pulpectomy in primary and permanent teeth: a systematic review and meta-analysis', *Journal of Clinical and Experimental Dentistry*, 16(9), pp. 1120-1128. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC11470449/>
- 2 Kaul, S., Kumar, A., Jasrotia, A. et al. (2021) 'Comparative Analysis of Biodentine™, Calcium Hydroxide, and 2% Chlorhexidine with Resin-modified Glass Ionomer Cement as Indirect Pulp Capping Materials in Young Permanent Molars', *Journal of Contemporary Dental Practice*, 22(5), pp. 511–516. Available at: <https://pubmed.ncbi.nlm.nih.gov/34318769/>
- 3 Internal market research data.
- 4 Duncan HF, Galler KM, Tomson PL, et al. European Society of Endodontology position statement: Management of deep caries and the exposed pulp. *Int Endod J*. 2019;52(7):923–934.. Available at: <https://onlinelibrary.wiley.com/doi/epdf/10.1111/iej.13080>
- 5 American Association of Endodontists (2021) *Vital Pulp Therapy: AAE Position Statement*. Available at: https://www.aae.org/wp-content/uploads/2021/05/VitalPulpTherapyPositionStatement_v2.pdf
- 6 Watson, T.F., Atmeh, A.R., Sajini, S. et al. (2014) 'Present and future of glass-ionomers and calcium-silicate cements as bioactive materials in dentistry: biophotonics-based interfacial analyses in health and disease', *Dental Materials*, 30(1), pp. 50–61. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC3885799/>
- 7 Camilleri, J. (2011) 'Characterization and hydration kinetics of tricalcium silicate cement for use as a dental biomaterial', *Dental Materials*, 27(8), pp. 836–844. Available at: <https://www.sciencedirect.com/science/article/abs/pii/S0109564111001151>
- 8 Camilleri, J. (n.d.) *Biodentine™: The dentine in a capsule or more*. Septodont Learning. Available at: <https://www.septodontlearning.co.uk/sites/default/files/Biodentine-The-dentine-in-a-capsule-or-more-by-Josette-Camilleri-compressed.pdf>
- 9 Internal data: Biodentine™ Scientific file (2011). Pages 20-23. Pages 24-28.
- 10 Hashem, D., Mannocci, F., Patel, S. et al. (2015) 'Clinical and radiographic assessment of the efficacy of calcium silicate indirect pulp capping: a randomized

controlled clinical trial', *Journal of Dental Research*, 94(4), pp 562-568. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC4485218/>

- 11 Nowicka, A., Lipski, M., Parafiniuk, M. *et al.* (2013) 'Response of human dental pulp capped with Biodentine and mineral trioxide aggregate', *Journal of Endodontics*, 39(6), pp. 743–747. Available at: [https://www.jendodon.com/article/S0099-2399\(13\)00133-7/abstract](https://www.jendodon.com/article/S0099-2399(13)00133-7/abstract)
- 12 Kaur, M., Singh, H., Dhillon, J.S. *et al.* (2017) 'MTA versus Biodentine™: Review of literature with a comparative analysis', *Journal of Clinical and Diagnostic Research*, 11(8), pp. ZG01–ZG05. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC5620936/>
- 13 Sajini, S., Atmeh, A., Banerjee, A. *et al.* (2022) 'Glass-ionomer and calcium silicate-based cements interactions with human dentine in health and disease: Two-photon fluorescence microscopy and Raman spectroscopy analysis', *Dental Materials*, 38(11), pp. 1710–1720. Available at: <https://www.sciencedirect.com/science/article/abs/pii/S0109564122002536?via%3Dihub>
- 14 Atmeh, A.R., Chong, E.Z., Richard, G. *et al.* (2012) 'Dentin-cement interfacial interaction: calcium silicates and polyalkenoates', *Journal of Dentistry*, 91(5), pp. 454–459. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC4077527/>
- 15 Elbanna, A., Atta, D., Sherief, D.I. (2022) 'In vitro bioactivity of newly introduced dual-cured resin-modified calcium silicate cement', *Dental Research Journal (Isfahan)*, 28(19), p. 1. Available at: <https://pubmed.ncbi.nlm.nih.gov/35308449/>
- 16 Nowicka, A., Wilk, G., Lipski, M. *et al.* (2015) 'Tomographic evaluation of reparative dentin bridge formation after direct pulp capping with calcium hydroxide, MTA, Biodentine, and dentin bonding system in human teeth', *Journal of Endodontics*, 41(8), pp. 1234–1240. Available at: <https://pubmed.ncbi.nlm.nih.gov/26031301/>
- 17 Gong, V., Franca, R. (2017) 'Nanoscale chemical surface characterization of four different types of pulp-capping materials', *Journal of Dentistry*, 58, pp. 11-18. Available at: <https://pubmed.ncbi.nlm.nih.gov/28043847/>
- 18 About, I. (2022) *Biodentine: Properties and Clinical Applications*. Switzerland: Springer. URL: <https://link.springer.com/book/10.1007/978-3-030-80932-4>
- 19 Lundy, F.T., Giraud, T., El-Karim, I.A. *et al.* (2022) 'Biodentine™ in inflammation and pain control', in About, I. (ed.) *Biodentine: Properties and Clinical Applications*. Switzerland: Springer, pp. 51–66. URL: <https://link.springer.com/book/10.1007/978-3-030-80932-4>

- 20 Niranjan, B., Shashikiran, N.D., Singla, S. *et al.* (2016) 'A comparative microleakage evaluation of three different base materials in Class I cavity in deciduous molars in sandwich technique using dye penetration and dentin surface interface by scanning electron microscope', *Journal of Indian Society of Pedodontics and Preventive Dentistry*, 34(4), pp. 324–330. Available at: <https://pubmed.ncbi.nlm.nih.gov/27681395/>
- 21 Abdelmegid, F. Y., Salama, F. S., Al-Mutairi, W. M., *et al.* (2017) 'Effect of different intermediary bases on microleakage of a nanocomposite interface in Class II box cavities in primary teeth: An in vitro study', *International Journal of Artificial Organs*, 40(2), pp. 82-87. Available at: <https://pubmed.ncbi.nlm.nih.gov/28222209/>